

A Donaldson Company

A WORLD LEADER IN FUME **EXTRACTION TECHNOLOGY**

DentalPRO Xtract 200





Dust and fume extraction cabinet for dental laboratory applications.

BOFA's DentalPRO Xtract 200 dust and fume extraction cabinet is designed for the removal of dust and fumes generated during dental hand finishing or monomer mixing applications.

This unit is ideally suited to manual benchtop applications and provides technicians with the flexibility to work without restriction in the knowledge that dust particles and fumes are being removed to maintain a safe working environment.

This system combines maximum technician protection in a small footprint, with scalable and flexible options.

Technology



ProTECT service plan



quality standard

Key features of the DentalPRO Xtract 200

Easily replaced screen protector ensures clear visibility at all times

Standard

Side entry grommets for mains powered hand finishing units Standard

Strong rigid construction Standard

DentalPRO 250 filtration system Optional

DentalPRO Mono filtration system Optional

Internal LED strip light

Standard

Removable collection tray

Standard

Powder coated white

Standard


DentalPRO 400 filtration system

Optional

Contact BOFA at www.donaldsonbofa.com/contact/



Technical specification



Airflow through filters



Contaminated air



LED lighting



Technical data					
	EU	US			
Dimensions	355 x 540 x 400mm	13.98 x 21.26 x 15.75"			
Opening dimensions (HxW)	315 x 540mm	12.4 x 21"			
Recommended face velocity	0.5m/s	150 f/min			
Approvals	UKCA and CE	UKCA and CE			

Cabinet part numbers		
Model	Part number	
DentalPRO Xtract 200	D0946A0000	

Hose kit part numbers (needed if ordering extraction system)		
Description	Part number	
Hose kit to connect extraction systetem to DentalPRO Xtract cabinet	A1020068	

Extraction system part numbers					
DentalPRO 400					
Model	Option	Voltage	Part number		
DentalPRO 400 powder coated unit twin port 75mm lid		90V - 257V	DL0944A0000		
DentalPRO 400 powder coated unit twin port 75mm lid	24V stop / start override silter change signal	90V - 257V	DL0944A0001		
Unit options					
Foot switch (on / off operation)			A1070184		
DentalPRO 250					
Model		Voltage	Part number		
DentalPRO 250 powder coated unit single port 75mm lid		230V	DL0142A0000		
DentalPRO 250 powder coated unit single port 75mm lid		115V	DL0141A0000		
DentalPRO Mono					
		w h	D		
Model		Voltage	Part number		
DentalPRO Mono powder coated unit single port 75mm lid		230V	DL2043A0000		
DentalPRO Mono powder coated unit single port 75mm lid		115V	DL2041A0000		

Other hose kits and filters are available.

Datasheet correct at time of publishing.

Where applicable, the carbon used in BOFA units is capable of removing a wide range of VOCs, however it is the responsibility of the user to ensure the carbon is suitable for their application. For specific applications, please contact us for details.

Important Notice: Many factors beyond the control of BOFA can affect the use and performance of BOFA products in a particular application, including the conditions under which the product is used. Since these factors are uniquely within the user's knowledge and control, it is essential the user evaluate the products to determine whether the product is fit for the particular purpose and suitable for the user's application. All products, product specifications, availability and data are subject to change without notice, and may vary by region or country.

Think before you print! Please consider the environment before printing this document.